

<b>PB2 Water Quality Monitor Adjustment to Alarm on %-Rejection for Nuisance Alarms</b>		TSB 2012003
Created/Last Revised Date 04/27/12	Last Reviewed Date 04/27/12	Page 1 of 3

**TECHNICAL SERVICE BULLETIN**

## ISSUE

On occasion feed water quality changes for various reasons. This is sometimes due municipalities changing water sources, modifying pretreatment, or construction and repair to main water lines. These can cause erratic feed TDS changes and fluctuating pH.

The PB2 is set at the factory to alarm on %-Rejection. In this setting the water quality is more prone to a nuisance alarm due to the above mentioned changes in feed water quality. The PB2 can alternately be set to allow the alarm set-point to be set on Product TDS instead of %-Rejection.

**WARNING**

- 1. ELECTRICAL HAZARD: Be careful where and what is touched when making the changes described below inside the internal control box.**
- 2. The following changes should not be done while a hemodialysis treatment is being performed.**

## SOLUTION

### PRE-CHANGE VERIFICATION


1. Verify the feed and product TDS values displayed on the Water Quality Monitor with a hand-held TDS meter. - If out of calibration see the Operators Manual for the calibration procedure. This manual is available online and can be downloaded from [www.betterwater.com/support](http://www.betterwater.com/support).

### DIPSWITCH CHANGE

1. Turn the device **OFF**.
2. Unplug from the electrical receptacle
3. Open the **External Top Panel** and the remove the cover from the internal **Control Box** and locate the control board containing the...

... block of two **Dipswitches**...

...the **CAL-UP** and **CAL-DOWN** buttons



View from the BACK of the PB2



<b>PB2 Water Quality Monitor Adjustment to Alarm on %-Rejection for Nuisance Alarms</b>		TSB 2012003
Created/Last Revised Date 04/27/12	Last Reviewed Date 04/27/12	Page 3 of 3

8. Replace the cover on the internal **Control Box**, and close the outer **External Top Panel** cover.

**NOTE:** This same product set-point procedure applies if the Water Quality Monitor's Right Dipswitch is set to **OFF:µS** which will display in micro Siemens rather than parts-per-million if the switch is set to **ON:PPM**.