









Dual Tank Bicarb – Short Profile









RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
Propeller Shaft - 316 Stainless Steel	EQBICB00471	
Propeller - Polypropylene	EQBICB01923	
Propeller Set Screw - Stainless Steel; one per propeller	EQBICB01920	
Shaft Coupling - Stainless Steel	EQBICB00472	
Coupling Set Screw - Stainless Steel; four per coupling	HWSCSS01920	
Mixer Motor - 1/4 HP, 115 VAC	EQBICB01868	
Mixer-Motor Power Cord - for Mix/Delivery Bicarb units	ELPCOO00200	
MD 100RLT Pump - 1/3 HP, 1 Phase, 115 VAC	EQPUIW00453	

Pictures do not reflect the size of the item in relation to the other pictures








Dual Tank Bicarb – Short Profile

RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
1000 FT Loop Pump - 3/4 HP, 1 Phase, 115 VAC	EQPUGR00600	
High Level Proximity Sensor Subassembly, Tank 1	EQSUBTNK1HL	
Low Level Proximity Sensor Subassembly, Tank 1	EQSUBTNK1LL	
High Level Proximity Sensor Subassembly, Tank 2	EQSUBTNK2HL	
Low Level Proximity Sensor Subassembly, Tank 2	EQSUBTNK2LL	
Flow Control, 2.0 GPM PVC Sch-80	PLFCS802001	
Check Valve, 1/2", 316 Stainless-Steel	PLVASS00846	
Flow-Switch, 3/4", PVC-Clear	PLFSS800444	

Pictures do not reflect the size of the item in relation to the other pictures





Dual Tank Bicarb – Short Profile

DESCRIPTION	RELATED REPLACEMENT PARTS PART#	PICTURE
Volara Foam for Tank Lid, .25" TK x 1.5" W * 6 ft for 60 gallon Tanks * 10 ft for 100 gallon Tanks	EQBICB01600	
1/2" Red Handle Valve, Asahi * Valve B21	PLVAS800167-A	
1/2" Blue Handle Valve, Asahi * Valves B6, B7	PLVAS800198-A	
3/4" Red Handle Valve, Asahi * Valve B8	PLVAS800169-A	
3/4" Blue Handle Valve, Asahi * Valves B11, B12	PLVAS800170-A	
1" Red Handle Valve, Asahi * Valves B13, B14	PLVAS800174-A	
1/2" Red Handle, Compact Valve * Valves B4, B5	PLVAS8C0004	

Pictures do not reflect the size of the item in relation to the other pictures

Dual Tank Bicarb – Short Profile

RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
Bicarb Unit Control Box	EQASSYBCB01854	
Fuse, 3 Amp, Little-Fuse <i>* For models manufactured before mid-September 2013</i>	ELLFFS00834	
Fuse, 2 Amp, Little-Fuse <i>* For models manufactured after mid-September 2013</i>	ELLFFS00832	
Bicarb Remote Alarm Box	EQASSYBCB01709	
Concentrate Regulator Calibration Kit <i>* For calibrating regulators in Floor-Valve-Boxes and Panels</i> <i>* Blue wand is for bicarbonate; Red wand is for acid</i>	EQASSYCC01	

Pictures do not reflect the size of the item in relation to the other pictures