



Better Water LLC

## Dual Tank Bicarb-Short Profile, Central Mix and Delivery

### Dual Tank Bicarb-Short Profile, Central Mix & Delivery



60 Gallon Bicarb

- Better Water LLC has two different sized models: 60 Gallon Tank and 100 Gallon Tank
- Provide an effective means to mix and distribute bicarbonate solution for use in hemodialysis applications.
- Can mix and deliver from either tank.
- Simple operation via a control box.
- Short Profile (*only 54" from floor to top of tank*) for easier access.

### Standard Features

- \* Rugged, Long-Lasting, Polyethylene Tanks
- \* Tanks have graduated gallon marks
- \* Hinged Lid
- \* Mixer-Motors
- \* Front Access Panel
- \* Distribution Pump
- \* Built-in header for loop distribution
- \* Control Box
- \* Jug Fills and Sample Ports for Each Tank
- \* Rinse Down Hose Connects to Spray-Balls in each Tank



100 Gallon Bicarb

### SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. May 2016



Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.






Better Water LLC

# Dual Tank Bicarb, Central Mix and Delivery

## Models and Specifications

Model	EQASSY60BC-2TNK-SP	EQASSY100BC-2TNK-SP
Total Capacity	60 Gallons	100 Gallons
Dimensions	Height 70" x Length 70" x Width 33" 54" (floor to top of tank)	Height 70" x Length 81" x Width 40" 54" (floor to top of tank)
Operating Weight	1350 lbs	2100 lbs
Tank Material	Polyethylene	
Electrical Requirements	110 VAC, 20 amps, 60 Hz, GFCI Protected	
Control Voltage	24 VAC	
Loop Connection	3/4" PVC Sch-80 Female Socket, solvent weld	
RO Water Connection	1" Stainless-Steel, Female Threaded Socket	
Drain Connection	1 1/2" PVC Sch-80 Female Socket, solvent weld	



## SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. May 2016



Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.

