

# BETTER WATER LLC

## KB Filter Rack

### KB FILTER RACKS



Better Water LLC KB filter racks have .03 $\mu$  pleated Filter Elements that are enclosed in filter housings to process RO Product water or DI polished RO product water. This filter assembly is not intended for any other application.

The .03 $\mu$  Filter is a pyrogen filter and will reject (capture) bacteria, pyrogens, and endotoxins.

The .03 $\mu$  filter membranes in this assembly are of the highest quality known.

Each .03 $\mu$  filter assembly is built to the specifications of the application.

The application will determine the number of filters on each assembly.

### Standard Features

- .03  $\mu$  membranes
- Pressure gauges
- 1,2,3 and 4 filters are available
- Sample port
- Rugged steel frame
- Powder coated

### Solutions. Service. Support.

Better Water LLC 698 Swan Dr Smyrna, Tn 37167 (615) 355-6063 and at [www.betterwater.com](http://www.betterwater.com)

Equipment changes are at the sole discretion of BW,LLC. Brochure photographs & specs may not reflect recent changes



# KB Filter Rack



Models and Specifications	
<b>Model (example only)</b>	<b>Eqassy0.75kB3.03</b>
<b>Capacity GPM</b>	12
<b># of membranes</b>	3
<b>Membranes</b> (1,2,3 and 4 membranes assemblies available)	.03 pleated filters
<b>Overall dimensions</b>	Typically 24" w x 22"d x 36" h wall mounted
<b>Operating weight</b>	40 lbs
<b>Incoming water Requirements</b>	R.O. or D.I polished water
<b>Gauges</b>	Stainless steel accurate to within 1%,
<b>Disinfecting</b>	Compatible with 1% Peracetic acid (Minnicare) and bleach (500ppm)

***BW LLC carries over 2000 parts, accessories and consumables in stock, and ready to ship including:***

- Carbon Filters (provides 10 Minute Empty Bed Contact Time)
- Duplex & Simplex Softener
- Storage Tanks
- Brass pressure regulator
- Blend valve
- City boost pump
- Thermometer
- Floor valve boxes
- Remote alarms
- & More!



Central Water Systems	Central Delivery Systems	Reverse Osmosis Units	Portable RO Units
Pre-Treatment Racks	Post Treatment Racks	D.I. Exchange Accessories	Consumables
<b>Hemodialysis</b>	<b>Medical</b>	<b>Pharmaceutical</b>	<b>Laboratories</b>
		<b>Commercial</b>	<b>Industrial</b>