

### Post Treatment Rack



Better Water, Inc has the solution to all of your post Reverse Osmosis and final filtration needs. The Post Treatment Frame utilizes a single frame to accommodate re-pressurization pumps, D.I. connections, D.I. resin trap, final filtration, and all the necessary controls and safety features. Each unit has the same features and safety considerations of Better Water's standard post treatment equipment. Instead of being manufactured on a job site or in an operational clinic, they can now be manufactured in a controlled environment. This allows for unparalleled quality and consistency during production of what is one of the most important parts of the water treatment system - *Final Filtration*.

The units come ready for easy installation. Simply plumb to the storage tank, the R.O. Water loop feed, and the drain, then it's ready to go. This makes them the perfect fit for the new water treatment system installations or retrofitting to an existing facility. Options include Resistivity or Conductivity Monitors & Dual Re-pressurization Pumps.

### Standard Features

- Re-pressurization pump(s)
- D.I. connections
- D.I. resin trap
- Final filtration
- Flow switch
- Sample port
- Gauges
- Simple operation
- Rugged metal frame construction
- Powder coated finish



**Solutions. Service. Support**

**[www.betterwater.com](http://www.betterwater.com)**

Equipment changes are at the sole discretion of BWI. Brochure photographs & specs may not reflect recent changes.



# POST TREATMENT

## MODELS AND SPECIFICATIONS

Part Number	EQPST 170-00-1AK	EQPST 305-00-1AK	EQPST 255-00-2AK	EQPST 340-00-2CK
Capacity GPM	10	12	15	20
# of Filters/Membranes	2	3	3	4
Filters/Membranes	4x40 spiral wound	20" hollow fiber	4X40 spiral wound	4x40 spiral wound
Operating Weight	320lbs	290lbs	425lbs	510lbs
Overall Dimensions	25" W x 45" D x 72" H			
Power Requirements	208 vac/3phase/60Hz (Single Phase available)			
Number of Pumps	1	1	2	2
Incoming water Requirements	Reverse Osmosis Water			
Gauges	All Stainless Steel construction, 1.5 % accurate, located pre and post the unit			
Other standard features	Sample port and gauges, resin trap, flow switch, compatible with 1 % Peracetic acid and bleach for disinfecting			

*Custom configurations available upon request*

**BWI carries over 2000 parts, accessories and consumables in stock, and ready to ship including:**

Carbon Filters (provides 10 Minute Empty Bed Contact Time)  
 Storage Tanks  
 Duplex & Simple Softener  
 Brass pressure regulator  
 Blend valve  
 City boost pump  
 Thermometer  
 Floor valve boxes  
 Remote alarms



Central Water Systems    Central Delivery Systems    Reverse Osmosis Units    Portable RO Units

Pre-Treatment Racks    Post Treatment Racks    D.I. Exchange    Accessories    Consumables

Hemodialysis    Medical    Pharmaceutical    Laboratories    Commercial    Industrial

Better Water, Inc  
 698 Swan Drive Smyrna, TN 37167  
 Phone: 615.355.6063    Fax: 615.355.6065  
 www.betterwater.com