




Better Water LLC

# Single Tank Bicarb, Central Mix and Delivery



Wall-Mounted Header



Control Box

## Standard Features

- \* Rugged, Long-Lasting, Polyethylene Tank
- \* Tanks have graduated gallon marks
- \* Hinged Lid
- \* Mixer-Motor
- \* Jug-Fill Valve
- \* Distribution Pump
- \* Wall-Mounted Header for Loop Connectivity
- \* Control Box

## Single Tank Bicarb, Central Mix & Delivery

- Better Water LLC has two different sized models based on your individual needs:
  - 60 Gallon Tank
  - 100 Gallon Tank
- Both provide an effective means to mix and distribute bicarbonate solution for use in hemodialysis applications.
- A wall-mounted Header provides the connection point to the distribution loop.
- Simple operation via a wall-mounted control box.



60 Gallon Bicarb



100 Gallon Bicarb

## SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. Mar 2014



*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*






Better Water LLC

# SingleTank Bicarb, Central Mix and Delivery

## Models and Specifications

Model	EQASSY60BC-1TNK	EQASSY100BC-1TNK
Total Capacity	60 Gallons	100 Gallons
Tank Material	Polyethylene	Polyethylene
Electrical Requirements	110 VAC, 20 amps, 60 Hz, GFCI Protected	110 VAC 20 amps, 60 Hz, GFCI Protected
Dimensions	Height 6' x Length 2' 6" x Width 3' 6"	Height 6' 2" x Length 3' 3" x Width 4'
Operating Weight	790 lbs	1220 lbs
Photo		

## SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. Mar 2014



*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*