



Better Water LLC

## Single Tank Bicarb, Mix-Only

### Single Tank Bicarb, Mix-Only



100 Gallon Bicarb

- Better Water LLC has two different sized models based on your individual needs:
  - 60 Gallon Tank
  - 100 Gallon Tank
- Both provide an effective means to mix bicarbonate solution for use in hemodialysis applications.
- Easy access to the tank is provided via a hinged lid on top.
- Simple mixer operation via an ON/OFF Switch.

### Standard Features

- \* Rugged, Long-Lasting, Polyethylene Tank
- \* Tanks have graduated gallon marks
- \* Hinged Lid
- \* Check Valve on Fill Port
- \* Mixer-Motor
- \* Jug-Fill Valve
- \* Drain Valve



60 Gallon Bicarb

## SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. Mar 2014



*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*







Better Water LLC

# SingleTank Bicarb, Mix-Only

## Models and Specifications

Model	EQASSY60BICARB	EQASSY100BICARB
Total Capacity	60 Gallons	100 Gallons
Tank Material	Polyethylene	Polyethylene
Electrical Requirements	110 VAC, 20 amps, 60 Hz, GFCI Protected	110 VAC 20 amps, 60 Hz, GFCI Protected
Dimensions	Height 6' x Length 2' 6" x Width 3' 6"	Height 6' 2" x Length 3' 3" x Width 4'
Operating Weight	790 lbs	1220 lbs
Photo		

## SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. Mar 2014



*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*