



Better Water LLC

Cartridge Based Pyrogen Filter System (2-4 filters)



System in picture shown
with 3 filter cartridges

- Better Water LLC's Cartridge Based Filter System utilizes filter cartridges enclosed in housings to process RO/DI water and is wall-mounted. This filter assembly is not intended for any other applications.
- The Pyrogen filter will reject (*capture*) bacteria, pyrogens, and endotoxins.
- There are two filter options; .03 micron and .05 micron
- Available with 2, 3, and 4 filter cartridges. The specific application will determine the number of filters required.
- Includes inlet and outlet pressure gauges, drain port, and sanitary sample ports.

Standard Features

- * .03 micron or .05 micron filter cartridge options
- * 2, 3, and 4 filter cartridge options
- * Pressure Gauges (*inlet and outlet*)
- * Sanitary Sample Ports (*inlet and outlet*)
- * Rugged steel frame
- * Powder-coated

SOLUTIONS * SERVICE * SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065
Website www.betterwater.com

rev. Apr 2017



Equipment changes are at the sole discretion of Better Water LLC.
Brochure photographs and specifications may not reflect recent changes.






Better Water LLC

Cartridge Based Pyrogen Filter System (2-4 filters)

Models and Specifications			
.03 micron Filter Models	EQASSY0.75KB2.03	EQASSY0.75KB3.03	EQASSY0.75KB4.03
.05 micron Filter Models	EQASSY0.75KB2.05	EQASSY0.75KB3.05	EQASSY0.75KB4.05
Number of Filters	2	3	4
Capacity GPM	.03 filter = 20 gpm .05 filter = 12 gpm	.03 filter = 30 gpm .05 filter = 18 gpm	.03 filter = 40 gpm .05 filter = 24 gpm
Pyrogen Filter	.03 micron, pleated filter, 2.5" x 20" .05 micron, hollow fiber filter, 2.5" x 20"		
Overall Dimensions	24" Depth x 36" Width x 36" Height <i>(max for 4 membranes)</i>		
Operating Weight	20 lbs	40 lbs	60 lbs
Incoming Water Requirements	RO or DI water		



SOLUTIONS * SERVICE * SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065
Website www.betterwater.com

rev. Apr 2017



Equipment changes are at the sole discretion of Better Water LLC.
Brochure photographs and specifications may not reflect recent changes.

