



Better Water LLC

## Dual Tank Bicarb, Central Mix and Delivery

### Dual Tank Bicarb, Central Mix & Delivery



60 Gallon Bicarb

- Better Water LLC has two different sized models based on your individual needs:
  - 60 Gallon Tank
  - 100 Gallon Tank
- Both provide an effective means to mix and distribute bicarbonate solution for use in hemodialysis applications.
- Can mix and deliver from either or both tanks.
- Simple operation via a control box.

### Standard Features

- \* Rugged, Long-Lasting, Polyethylene Tank
- \* Tanks have graduated gallon marks
- \* Hinged Lid
- \* Mixer-Motor
- \* Jug-Fill Valve
- \* Distribution Pump
- \* Built-in header for loop distribution
- \* Control Box



100 Gallon Bicarb

### SOLUTIONS \* SERVICE \* SUPPORT

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)

rev. Oct 2014





*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*





Better Water LLC

# Dual Tank Bicarb, Central Mix and Delivery

Models and Specifications		
Model	EQASSY60BC-2TNK	EQASSY100BC-2TNK
Total Capacity	60 Gallons	100 Gallons
Tank Material	Polyethylene	Polyethylene
Electrical Requirements	110 VAC, 20 amps, 60 Hz, GFCI Protected	110 VAC 20 amps, 60 Hz, GFCI Protected
Control Voltage	24 VAC	24 VAC
Dimensions	Height 78" x Length 70" x Width 30" 60 1/2" (floor to top of tank)	Height 80" x Length 81" x Width 37" 63" (floor to top of tank)
Operating Weight	1350 lbs	2100 lbs
Loop Connection	3/4" PVC Sch-80 Female Socket, solvent weld	3/4" PVC Sch-80 Female Socket, solvent weld
RO Water Connection	1" Stainless-Steel, Female Threaded Socket	1" Stainless-Steel, Female Threaded Socket
Drain Connection	1 1/2" PVC Sch-80 Female Socket, solvent weld	1 1/2" PVC Sch-80 Female Socket, solvent weld
Photo		

## SOLUTIONS \* SERVICE \* SUPPORT

rev. Oct 2014

Better Water LLC; 698 Swan Dr; Smyrna, TN 37167; Phone (615) 355-6063, Fax (615) 355-6065  
Website [www.betterwater.com](http://www.betterwater.com)



*Equipment changes are at the sole discretion of Better Water LLC.  
Brochure photographs and specifications may not reflect recent changes.*