









Dual Tank Bicarb-Standard










RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
Propeller Shaft - 316 Stainless Steel	EQBICB00471	
Propeller - Polypropylene	EQBICB01923	
Propeller Set Screw - Stainless Steel	EQBICB01920	
Shaft Coupling - Stainless Steel	EQBICB00472	
Coupling Set Screw - Stainless Steel	HWSCSS01920	
Mixer Motor - 1/4 HP, 115 VAC	EQBICB01868	
Mixer-Motor Power Cord - for Mix/Delivery Bicarb units	ELPCOO00200	
MD 100RLT Pump - 1/3 HP, 1 Phase, 115 VAC	EQPUIW00453	

Pictures do not reflect the size of the item in relation to the other pictures

Dual Tank Bicarb-Standard

RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
1000 FT Loop Pump - 3/4 HP, 1 Phase, 115 VAC	EQPUGR00600	
High Level Proximity Sensor Subassembly, Tank 1 <i>* For models manufactured from February 2010</i>	EQSUBTNK1HL	
Low Level Proximity Sensor Subassembly, Tank 1 <i>* For models manufactured from February 2010</i>	EQSUBTNK1LL	
High Level Proximity Sensor Subassembly, Tank 2 <i>* For models manufactured from February 2010</i>	EQSUBTNK2HL	
Low Level Proximity Sensor Subassembly, Tank 2 <i>* For models manufactured from February 2010</i>	EQSUBTNK2LL	
High Level Float Switch, Tank 1 <i>* For models manufactured before February 2010</i>	EQSUBBICBTNK1FL	
Low Level Float Switch, Tank 1 <i>* For models manufactured before February 2010</i>	EQSUBTNK1FL LOW	
High Level Float Switch, Tank 2 <i>* For models manufactured before February 2010</i>	EQSUBBICBTNK2FL	
Low Level Float Switch, Tank 2 <i>* For models manufactured before February 2010</i>	EQSUBTNK2FL LOW	

Pictures do not reflect the size of the item in relation to the other pictures

Dual Tank Bicarb-Standard

RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
Bicarb Unit Control Box	EQASSYBCB01854	
Bicarb Remote Alarm Box	EQASSYBCB01709	
Concentrate Regulator Calibration Kit <i>* For calibrating regulators in Floor-Valve-Boxes and Panels</i> <i>* Blue wand is for bicarbonate; Red wand is for acid</i>	EQASSYCC01	
Flow Control, 2.0 GPM PVC Sch-80	PLFCS802001	
Check Valve, 1/2", 316 Stainless-Steel	PLVASS00846	
Fuse, 3 Amp, Little-Fuse <i>* For models manufactured before mid-September 2013</i>	ELLFFS00834	
Fuse, 2 Amp, Little-Fuse <i>* For models manufactured from mid-September 2013</i>	ELLFFS00832	
Flow-Switch, 3/4", PVC-Clear	PLFSS800444	

Pictures do not reflect the size of the item in relation to the other pictures

Dual Tank Bicarb-Standard



RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
Volara Foam for Tank Lid, .25" TK x 1.5" W * 6 ft for 60 gallon Tanks * 10 ft for 100 gallon Tanks	EQBICB01600	
1/2" Red Handle Valve, Plast-O-Matic * For models manufactured before mid-July 2013 * Valve B21	PLVAS800167	
1/2" Red Handle Valve, Asahi * For models manufactured from mid-July 2013 * Valve B21	PLVAS800167-A	
3/4" Red Handle Valve, Plast-O-Matic * For models manufactured before mid-July 2013 * Valve B8	PLVAS800169	
3/4" Red Handle Valve, Asahi * For models manufactured from mid-July 2013 * Valve B8	PLVAS800169-A	
3/4" Blue Handle Valve, Plast-O-Matic * For models manufactured before mid-July 2013 * Valves B6, B7, B9, B10, B15	PLVAS800170	
3/4" Blue Handle Valve, Asahi * For models manufactured from mid-July 2013 * Valves B6, B7, B9, B10, B15	PLVAS800170-A	
1" Blue Handle Valve, Plast-O-Matic * For models manufactured before mid-July 2013 * Valves B11, B12	PLVAS800172	
1" Blue Handle Valve, Asahi * For models manufactured from mid-July 2013 * Valves B11, B12	PLVAS800172-A	

Pictures do not reflect the size of the item in relation to the other pictures

Dual Tank Bicarb-Standard

RELATED REPLACEMENT PARTS

DESCRIPTION	PART#	PICTURE
1" Red Handle Valve, Plast-O-Matic <i>* For models manufactured before mid-July 2013</i> <i>* Valves B13, B14</i>	PLVAS800174	
1" Red Handle Valve, Asahi <i>* For models manufactured from mid-July 2013</i> <i>* Valves B13, B14</i>	PLVAS800174-A	

Pictures do not reflect the size of the item in relation to the other pictures